

# SNAP *Leishmania* Test

## Test accuracy

### Sensitivity and specificity of the SNAP *Leishmania* Test

Analyte	Reference method (IFA <sup>1</sup> )	SNAP <i>Leishmania</i> Test result		Total	Sensitivity (95% CL)
		+	-		Specificity (95% CL)
<i>Leishmania infantum</i> <sup>1</sup>	+	83	5	88	98,8% (93,5%–100%)
	-	1	268	269	98,2% (95,8%–99,4%)

**Table 1.** SNAP *Leishmania* Test versus reference methods.<sup>1</sup>

#### Reference

<sup>1</sup> Data on file at IDEXX Laboratories, Inc. Westbrook, Maine USA.